

How to prepare Autopen®

1. Remove the cap and unscrew the body from the cartridge holder.



2. Check to ensure that the internal plunger is fully retracted. Rewind the internal plunger by turning the body until it is fully retracted.



3. Load the insulin cartridge into the cartridge holder, metal cap first and screw the pen back together.



4. How to prime and inject...

Attach the Pentip and remove the protective covers.

5. Dial two units on the dose selector.



6. With the pen pointing upwards push and hold the release button until the starting line on the dose selector returns to the arrow on the body. Repeat these steps until a drop of insulin appears at the needle tip. Then dial the required dose.



7. Insert the needle, push and hold down the release button until the dose selector returns to the start position, count to 10 then remove the needle from the skin.



8. How to change the insulin cartridge...

Remove the Pentip and discard safely, unscrew the body and remove the cartridge.

9. Rewind the internal plunger by turning the body until it is fully retracted. Return to step 3.



If the dose selector stops before the starting line returns to the arrow, this indicates that you have not received a full dose.

The dose selector will show the number of units that must be injected from a new cartridge to complete the dose. Return to step 2.



These instructions are intended as a guide only. Full instructions are supplied with Autopen and must be read in full before using the device. Autopen must not be used without first consulting your medical adviser.