

Technical Specifications

- Operating voltage: with battery DC 12 V with power supply 100÷240 V ~ 50/60 Hz - 14 VDC 4,28A
- Power consumption: max 50 W
- Size: 134 x 378 x h 263 mm
- Weight: 5 kg
- Bottle capacity: 1000 or 2000 ml
- Flow: 30 air litres/min
- Sound level: 65 dBA at 1 m
- Operating cycles: continuous (45 minutes battery autonomy)
- MAX vacuum: 0 ± 85 kPa
- Internal battery: 12 VDC 4A h lead hermetic
- Battery autonomy: about 45 minutes, recharging time approx 6 hours
- Norms:
 - MDD 93/42/EEC (Class IIa - annex IX)
 - EN60601-1(IEC 60601-1) / EN60601-1-2 (IEC 60601-1-2)
 - UNI EN ISO 10079-1/ UNI CEI EN ISO 14971
 - Degree of protections against ingress of liquids: IP34d (IEC 529)
 - Class II- Type B