


Welch Allyn Green Series™ 900 Procedure Light



6 LEDs



Touchless On/Off
with Dimming Control



Disposable Sheath
Available—Part #52630

Whiter, Brighter, and Greener

The GS 900 Procedure Light features:

- Six 5,500° K LEDs that provide a whiter, brighter light than halogen
- Less energy consumption for reduced cost of ownership
- No bulb replacement for cost savings
- Combines outstanding durability, extraordinary maneuverability, unparalleled spot quality and multiple mounting options
- Dimming feature
- Minimal heat dissipation and high color temperature
- 5-year limited warranty



WelchAllyn®
Advancing Frontline Care™



Green Series Medical Exam Lights

Welch Allyn Green Series™ 900 Procedure Light

ORDERING INFORMATION

44900	GS 900 Procedure Light with Mobile Stand
44900-C	GS 900 Procedure Light with Ceiling Mount
44900-W	GS 900 Procedure Light with Wall Mount
52630	Disposable Sheath for GS 600/GS 900 (50/box, 5 boxes/case)

SPECIFICATIONS

# of LEDs	6 White LEDs
Color Temperature	5,500° K
Light Intensity @ Typical Working Distance	80,000 lux (7,432 footcandles)
LED Life	50,000 hours
Working Distance	24 in (60.9 cm) to 40 in (100 cm)
Spot Size at Working Distance	9 in (23 cm) to 11 in (28 cm)
Electrical	Input: 100-240 V~, 50 Hz-60 Hz, 1.00 A max
Power Cord	#18 AWG, 3-wire grounded cord set, 8' (2.5 m) long
Physical Dimensions	Weight: Total Wall Mount: 17.90 lbs (8.12 kg) Total Mobile Stand: 44.85 lbs (20.34 kg) Total Ceiling Mount: 24.50 lbs (11.11 kg) Luminaire: Height: 13.74 in (34.9 cm) Width: 18.03 in (45.8 cm) Depth: 4.02 in (10.2 cm) Wall to Luminaire Center: 68.9 in (175 cm) Wall Mount Max. Up/Down: 37.4 in (95 cm) Mobile Stand: Height: 68.9 in (175 cm) Mobile Stand Base: Length: 21.97 in (55.8 cm) Width: 22.24 in (56.5 cm) Height: 6.5 in (16.5 cm) Mobile Stand Max. Up/Down: 37.4 in (95 cm) Mobile Stand to Luminaire Center: 37.4 in (95 cm) Ceiling to Luminaire Center: 68.9 in (175 cm) Ceiling Mount Max. Up/Down: 37.4 in (95 cm)
Mounting Options	Wall Mount Mobile Stand Ceiling Mount
Control Box	Touchless On/Off Control with Dimming
Accessory Options	Disposable Sheath #52630
Warranty	5-year limited



IDEAL FOR USE IN THE FOLLOWING ENVIRONMENTS:

- AMBULATORY CARE CENTERS
- LABOR AND DELIVERY
- HOSPITALS—ER AND ICU
- OUTPATIENT FACILITIES
- OB/GYN
- URGENT CARE CENTERS

MOUNTING OPTIONS AVAILABLE:



To place an order, contact our Customer Service Department at 1.800.535.6663



WelchAllyn®
Advancing Frontline Care™

