

# Comparison Table

## Blood Pressure Monitors Fully Automatic Upper Arm

Which product is right for you? We made it simple for you to decide. The tables help you to review the many features of the OMRON products.



Model	MIT Elite Plus	MIT Elite	i-Q132	i-Q142	M10-IT	i-C10	M6 Comfort
Model Number	HEM-7301-ITKE7	HEM-7300-WE7	HEM-1010-E	HEM-1040-E	HEM-7080IT-E	HEM-7070-E	HEM-7223-E8
Intellisense™ Technology	●	●	●	●	●	●	●
AC-Adapter	E1600 (opt.) 3094298-6	E1600 (opt.) 3094298-6	E1000 (incl.) 3094100-9	E1600 (incl.) 3094298-6	S-Adpater (Opt.) 9515336-9	S-Adpater (Opt.) 9515336-9	S-Adpater (Opt.) 9515336-9
PC Connection with Software & Cable	●			●	●		
Supplied Cuff	Medium 22–32 cm	Medium 22–32 cm	Small/ Medium arm in style 17–32 cm	Medium/ Large arm in style 22–42 cm	Medium/ Large Preformed 22–42 cm	Medium/ Large Preformed 22–42 cm	Medium/ Large Preformed 22–42 cm
Optional Large Cuff Available	(MIT cuff)	(MIT cuff)		Built-in	Built-in	Built-in	Built-in
Optional Small Cuff Available			Built-in				
Memory Positions	90	90	84 p. user 2 users	84 p. user 2 users	84 p. user 2 users	84 p. user 2 users	100 p. user, 2 users
Three-fold Displays (Sys/Dia/Pulse)	● (with backlight)	●	●	●	●	●	●
Date/Time Indication	●	●	●	●	●	●	●
Large Display and Operation Buttons	●	●	●	●	●	●	●
Dual Check System							●
Cuff Wrapping Guide							●
Hypertension Indicator	Symbol	Symbol	Symbol	Symbol	Symbol	Symbol	Identification Bar
Irregular Pulse Detection	●	●	●	●	●	●	●
Averaging of 3 Values	●	●	●	●	●	●	●
Clinical Validation	IP	IP	Pending	IP	IP and BHS	IP and BHS	IP and BHS



M6W	705IT	705CP-II	M7	M3W	M3	M2
HEM-7213-E	HEM-759P-E	HEM-759P-E2	HEM-780-E	HEM-7202-E	HEM-7200-E2	HEM-7119-E
●	●	●	●	●	●	●
S-Adpater (Opt.) 9515336-9	S-Adpater (Opt.) 9515336-9	S-Adpater (Opt.) 9515336-9	S-Adpater (Opt.) 9515336-9	S-Adpater (Opt.) 9515336-9	S-Adpater (Opt.) 9515336-9	S-Adpater (Opt.) 9515336-9
	●	Optional				
Wide Range 22-42 cm	Medium 22-32 cm	Medium 22-32 cm	Medium/ Large Preformed 22-42 cm	Wide Range 22-42 cm	Medium 22-32 cm	Medium 22-32 cm
Built-in	●	●	Built-in	Built-in		
●	●	●				
100 p. user, 2 users	28	28	90	60 p. user, 2 users	60 p. user, 2 users	21
●	●	●	●	●	●	●
●	●	●	●	●	●	
●	●	●	●			
●				●	●	
●				●	●	●
Identification Bar				Identification Bar	Identification Bar	Identification Bar
●				●	●	●
●			●	●	●	
IP	IP and BHS	IP and BHS	IP and BHS	IP	IP	IP