

Class 4 Dual Multi Parameter Indicator

The Valisafe Class 4 Dual Multi Parameter indicators offer an economical yet effective way of monitoring your steam sterilization processes.

(i) With one accurate calibration world wide it is a secure and safe monitoring system yet simple and easy to use.

(ii) Unlike most other class 4 indicators available today the Valisafe indicator is non toxic - another sign of it's advanced technology.

- Conforms to EN ISO 11140–1 Class 4
- Class 4 Multi Parameter indicator - Steam
- Monitors all 121°C to 134°C steam sterilization cycles
- Clear accurate colour change
- Easy to interpret - easy to record
- For use in every tray or item
- Non-hazardous / Lead free according to EC Regulation 07/2006
- Strips dividable to make 500
- Laminated

Name:	Dual Multi Parameter Indicator
Classification:	EN-ISO 11140-1 (Class 4)
Operating temp:	121°C -134°C
Storage conditions:	10-30°C
RH:	30-70°C
Correct symbols:	In accordance with EN ISO 11140-1
Indicator:	Class 4
Disp/single use:	Yes
Coated:	Yes
Non-toxic/Lead free:	Yes, according to EC Regulation 1907/2006
Lot no.	Yes
Prod date/Exp date:	Yes
Operator:	Yes
Pass/fail result:	Yes
Colour change:	Blue to dark grey or black



Before
Original indicator



After
After total process



Order details:

V3502038

Dual Multi Parameter Indicator

250 Indicators / can be split (500)