

GIM-30724 Three LED light Sources Wall

Technical Specifications

Type of lighting: LED

Consumption (Watt): 8

Life of the bulbs (hours): > 60,000

Light intensity at 50 cm (lux): 60,000

Light intensity at 1 m (lux): 20,000

Colour temperature $\pm 5\%$ at 1 m (K): 4,500

Colour rendering index (Ra)%: 94

Light field size at 50 cm (\varnothing mm): 150

External \varnothing of the reflector (mm): 96

Colour: white

Norms: IEC 60601-1

Marks: CE 93/42

Weight (kg): 5

Operating voltage V 100-230

Operating Hz 50/60