



### TAKES YOUR MIND OFF THE PAIN

Pressing Unistik 3's 'Comfort Zone' of 8 raised dots against the skin sends a message of comfort to the brain. Then, when you click Unistik 3's release button, the lancet penetrates your skin, but the dominant Comfort Zone sensation means the pain 'message' to your brain is greatly reduced. Distracting the brain in this way makes sampling virtually pain-free.



Unistik 3 Comfort



Unistik 3 Normal



Unistik 3 Extra

## Ask your GP for Unistik 3 on prescription.

Three variants mean that Unistik 3 is suitable for all skin types.

Unistik 3 Variant	Number In Box	Product Code	PIP Code
Comfort	100	AT1042	317-2277
Comfort	200	AT1044	317-2285
Normal	100	AT1002	317-2301
Normal	200	AT1004	317-2293
Extra	100	AT1012	317-2269
Extra	200	AT1014	317-2251

For more information please call UK Customer Services on **01993 810052**, or visit [www.unistik3.com](http://www.unistik3.com)




[www.owenmumford.com](http://www.owenmumford.com)

UK & EXPORT  
Owen Mumford Ltd,  
Brook Hill, Woodstock,  
Oxford OX20 1TU, United Kingdom  
T : 01993 812021  
F : 01993 813466  
E : [info@owenmumford.co.uk](mailto:info@owenmumford.co.uk)

NORTH AMERICA  
Owen Mumford Inc,  
1755 West Oak Commons Court,  
Marietta, Georgia 30062, USA  
T : 800-421-6936  
F : 770-977-2866  
E : [info@owenmumfordinc.com](mailto:info@owenmumfordinc.com)

FRANCE  
Owen Mumford,  
BP 444, 27204 Vernon Cedex, France  
T : 02 32 51 88 70  
F : 02 32 21 96 95  
E : [contact@owenmumford-fr.com](mailto:contact@owenmumford-fr.com)

GERMANY  
Owen Mumford GmbH,  
Alte Häge 1,  
63762 Großostheim-Ringheim, Germany  
T : 06026/97750  
F : 06026/977529  
E : [info@owen-mumford.de](mailto:info@owen-mumford.de)



# How to make blood sampling less of a pain.

## Unistik® 3

Single-use safety lancets with Comfort Zone Technology®  
**Now available on prescription.**



AT1042/FOUN/0207/07

# Unistik 3<sup>®</sup>

Easy, comfortable blood sampling – **now available on prescription.**





With Unistik 3 there's now a convenient new way to test your blood glucose levels that's a lot less painful. That's because Unistik 3's unique Comfort Zone Technology<sup>®</sup> actually distracts your brain from feeling the prick of the needle.

## Safe and hygienic.

As Unistik 3 is a single-use safety lancet, you will reduce your risk of infection by using a fresh lancet every time. There's no need to load the lancet, and you never see it because it remains hidden inside the case. So there's also little chance of accidental injury before, or after, the test.


## There's a Unistik 3 safety lancet for everyone.

With three variants for different skin types, there's a Unistik 3 single-use safety lancet to suit every individual. And as each variant is a different colour, you'll always know that you're using the right one.



## Using Unistik 3 couldn't be simpler.

Simply twist off the protective sterile cap, press Unistik 3 on your finger, click the release button, test with any meter device and dispose. It's that easy.



For more information please call UK Customer Services on **01993 810052** or visit [www.unistik3.com](http://www.unistik3.com)