

Medema Limited. [3-9 Broomhill Road, London, SW18 4JQ. United Kingdom.]
 Telephone: 0044 (0) 208 8779258, Fax: 0044 (0) 208 8779888
 e-mail: sales@medema.co.uk, web: www.medema.co.uk

Resuscitator Breathing System for Single-Use [Resuscitator Breathing System for single-use]

Resuscitator Breathing System for single-use

- Disposable Mask #5
- NRV Complete with pop-off 60cm H2O for single use resuscitation
- Single use Adult 1500ml bag
- Intake valve complete
- Single use 1000ml reservoir
- Oxygen Tubing (Non-Klink) 2.0Mtr / 6mm OD, 4mm ID
- Ready-to-go design for any emergency applications
- Quality Waterproof Draw-string bag
- Well designed icons for Distinction
- Multinational languages for worldwise users



FULLY ASSEMBLED KIT
 Ready-to-go design for any emergency applications



COMPACT PACKING

[Product Code] Product Name	Image	Unit Price
<p>[SUGRB-0001] NEW Resuscitator Breathing System for single-use (Adult)</p>		<p>£10.00 (Ex. VAT)</p>
<p>[SUGRB-0002] NEW Resuscitator Breathing System for single-use (Child)</p>		<p>£9.50 (Ex. VAT)</p>