

The Unistik® 3 range

The Unistik 3 range of single use safety lancets has 4 models for use across all clinical environments:

Unistik® 3 *Comfort*

Gauge: 28G. Depth: 1.8mm

Average Blood Volume: 5 μ L-30 μ L.

Suggested Test Areas e.g. Ward blood glucose monitoring. HbA1c.

Suggested Patient Areas e.g. Diabetes Clinics. Paediatrics. Geriatrics. Home Use.

Box of 100: AT1042 / Box of 200: AT1044



Unistik® 3 *Normal*

Gauge: 23G. Depth: 1.8mm

Average Blood Volume: 30 μ L-75 μ L.

Suggested Test Areas e.g. Ward blood glucose monitoring. Cholesterol. HbA1c.

Suggested Patient Areas e.g. Adults. Geriatrics. Home Use. Blood transfusion centres.

Box of 100: AT1002 / Box of 200: AT1004



Unistik® 3 *Extra*

Gauge: 21G. Depth: 2.0mm

Average Blood Volume: 75 μ L-125 μ L.

Suggested Test Areas e.g. Cholesterol.

Haemoglobin (Hgb). Diagnostic kit tests such as cardiac enzymes.

Suggested Patient Areas e.g. Adults (tougher skin). Blood transfusion centres.

Box of 100: AT1012 / Box of 200: AT1014



Unistik® 3 *Neonatal & Laboratory*

Gauge: 18G. Depth: 1.8mm

Average Blood Volume: 200 μ L+.


Neonatal Suggested Test Areas e.g. Congenital Hypothyroidism. Phenylketonuria (PKU). Single heelstick tests such as Bilirubin, Glucose,

Haemoglobin. Sickle Cell disease. Other diseases.

Neonatal Suggested Patient Areas e.g.

Heelstick samples on neonates.

Box of 100: AT1052 / Box of 200: AT1054



Laboratory Suggested Test Areas e.g.

Blood Gases & Electrolytes. Chemistry Test Profiles. H.pylori.

Laboratory Suggested Patient Areas e.g.

ITU / ICU. CCU. Laboratories.